

# Overlay Report

Create a duplicate graph on top of an existing graph in a report template.



**Overlay\_Report** /PAGE\_NO= Select report page number or all pages /IDENTIFIER\_TYPE= Copy graph information for the identifiers types of either EVENT\_LABEL, TAG\_NAME or LINE\_COLOR /IDENTIFIER= Select the identifier name /NEW\_TAG= Select tag to overlay on current graph /NEW\_COLOR= Select graph color ! /OLD\_EVENT= ! /NEW\_EVENT= ; An example of using this command to compare data sets can be found in the [Comparing Data Sets Tutorial](#).

From:  
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:  
[https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:report\\_commands:overlay\\_report](https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:report_commands:overlay_report)

Last update: 2024/07/17 15:46

