

Create a duplicate graph on top of an existing graph in a report template.

[Overlay_Report.PNG](#)

Overlay_Report /PAGE_NO= Select report page number or all pages /IDENTIFIER_TYPE= Copy graph information for the identifiers types of either EVENT_LABEL, TAG_NAME or LINE_COLOR /IDENTIFIER= Select the identifier name /NEW_TAG= Select tag to overlay on current graph /NEW_COLOR= Select graph color ! /OLD_EVENT= ! /NEW_EVENT= ; An example of using this command to compare data sets can be found in the [Comparing Data Sets Tutorial](#).

From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:report_commands:overlay_report&rev=1718648554

Last update: 2024/06/17 18:22

