

create a duplicate graph on top of an existing graph in a report template.



overlay_report /page_no= select report page number or all pages /identifier_type= copy graph information for the identifiers types of either event_label, tag_name or line_color /identifier= select the identifier name /new_tag= select tag to overlay on current graph /new_color= select graph color ! /old_event= ! /new_event= ; an example of using this command to compare data sets can be found in the [comparing data sets tutorial](#).

From:

<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:

https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:report_commands:overlay_report&rev=1718801569

Last update: 2024/06/19 12:52

