

**Language:** English • français • italiano • português • español \*\*\*\*

==== Contents ====

- \* [1 Manually Set Report Options](#)
- \* [1.1 Page Layout](#)
- \* [1.2 Lab Information](#)
- \* [1.3 Page Layout](#)
- \* [1.4 Printing](#)

```
Set_Report_Options
! /PAGE_LAYOUT=LANDSCAPE
! /PAPER_SIZE=-
! /COLLECTOR=-
! /FACILITY_LINE_2=-
! /ADDRESS_LINE_1=-
! /ADDRESS_LINE_2=-
! /ADDRESS_LINE_3=-
! /PHONE=-
! /FAX=-
! /NOTES_LINE_1=-
! /NOTES_LINE_2=-
! /NOTES_LINE_3=-
! /LOGO=-
! /MARGIN_TOP=-
! /MARGIN_BOTTOM=-
! /MARGIN_LEFT=-
! /MARGIN_RIGHT=-
! /PAGE_NUMBER_POSITION=-
! /PAGE_NUMBER_ALIGNMENT=-
! /PAGE_TAG_LINE_2=-
! /PAGE_TAG_LINE_3=-
! /PAGE_TAG_LINE_4=-
;
```

## Manually Set Report Options

These parameters can be changed in the pipeline command, or changed manually by right clicking on the report page and selecting **Report Options**.

### Page Layout

#### Page Layout Tab:

**/PAGE\_LAYOUT=** Portrait or Landscape

**/PAPER\_SIZE=** Letter, Legal or A4

# Lab Information

The Lab Information can be added in the Laboratory Information tab.

## Laboratory Information Tab:

/COLLECTOR=  
/FACILITY\_LINE\_2=  
/ADDRESS\_LINE\_1=  
/ADDRESS\_LINE\_2=  
/ADDRESS\_LINE\_3=  
/PHONE=  
/FAX=  
/NOTES\_LINE\_1=  
/NOTES\_LINE\_2=  
/NOTES\_LINE\_3=  
/LOGO=

Note that to see the **lab information** in the report, the lab information must be “graphed”.  
By default facility line 1 is set to the laboratory name and this parameter cannot be modified.

# Page Layout

## Page Layout Tab:

1. Margins:  
/MARGIN\_TOP=  
/MARGIN\_BOTTOM=  
/MARGIN\_LEFT=  
/MARGIN\_RIGHT=
2. Page Number Format:  
/PAGE\_NUMBER\_POSITION=  
/PAGE\_NUMBER\_ALIGNMENT=

Note paper size and orientation can be set in the first two parameters of this pipeline command.

# Printing

## Printing Tab:

/PAGE\_TAG\_LINE\_2=  
/PAGE\_TAG\_LINE\_3=  
/PAGE\_TAG\_LINE\_4=

Retrieved from “”

From:  
<https://has-motion.com/wiki/> - **Wiki Documentation Site**

Permanent link:  
[https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:report\\_commands:set\\_report\\_options&rev=1718386185](https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:report_commands:set_report_options&rev=1718386185)

Last update: **2024/06/14 17:29**

