

Language: English • français • italiano • português • español ****

==== Contents ====

- * [1 Manually Set Report Options](#)
- * [1.1 Page Layout](#)
- * [1.2 Lab Information](#)
- * [1.3 Page Layout](#)
- * [1.4 Printing](#)

```
Set_Report_Options
! /PAGE_LAYOUT=LANDSCAPE
! /PAPER_SIZE=-
! /COLLECTOR=-
! /FACILITY_LINE_2=-
! /ADDRESS_LINE_1=-
! /ADDRESS_LINE_2=-
! /ADDRESS_LINE_3=-
! /PHONE=-
! /FAX=-
! /NOTES_LINE_1=-
! /NOTES_LINE_2=-
! /NOTES_LINE_3=-
! /LOGO=-
! /MARGIN_TOP=-
! /MARGIN_BOTTOM=-
! /MARGIN_LEFT=-
! /MARGIN_RIGHT=-
! /PAGE_NUMBER_POSITION=-
! /PAGE_NUMBER_ALIGNMENT=-
! /PAGE_TAG_LINE_2=-
! /PAGE_TAG_LINE_3=-
! /PAGE_TAG_LINE_4=-
;
```

Manually Set Report Options

These parameters can be changed in the pipeline command, or changed manually by right clicking on the report page and selecting **Report Options**.

Page Layout

Page Layout Tab:

/PAGE_LAYOUT= Portrait or Landscape

/PAPER_SIZE= Letter, Legal or A4

Lab Information

The Lab Information can be added in the Laboratory Information tab.

Laboratory Information Tab:

/COLLECTOR=
/FACILITY_LINE_2=
/ADDRESS_LINE_1=
/ADDRESS_LINE_2=
/ADDRESS_LINE_3=
/PHONE=
/FAX=
/NOTES_LINE_1=
/NOTES_LINE_2=
/NOTES_LINE_3=
/LOGO=

Note that to see the **lab information** in the report, the lab information must be “graphed”.
By default facility line 1 is set to the laboratory name and this parameter cannot be modified.

Page Layout

Page Layout Tab:

1. Margins:
/MARGIN_TOP=
/MARGIN_BOTTOM=
/MARGIN_LEFT=
/MARGIN_RIGHT=
2. Page Number Format:
/PAGE_NUMBER_POSITION=
/PAGE_NUMBER_ALIGNMENT=

Note paper size and orientation can be set in the first two parameters of this pipeline command.

Printing

Printing Tab:

/PAGE_TAG_LINE_2=
/PAGE_TAG_LINE_3=
/PAGE_TAG_LINE_4=

Retrieved from “”

From:

<https://has-motion.com/wiki/> - **Wiki Documentation Site**

Permanent link:

https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:report_commands:set_report_options&rev=1718644097

Last update: **2024/06/17 17:08**

