

```
Set_Report_Options
! /PAGE_LAYOUT=LANDSCAPE
! /PAPER_SIZE=-
! /COLLECTOR=-
! /FACILITY_LINE_2=-
! /ADDRESS_LINE_1=-
! /ADDRESS_LINE_2=-
! /ADDRESS_LINE_3=-
! /PHONE=-
! /FAX=-
! /NOTES_LINE_1=-
! /NOTES_LINE_2=-
! /NOTES_LINE_3=-
! /LOGO=-
! /MARGIN_TOP=-
! /MARGIN_BOTTOM=-
! /MARGIN_LEFT=-
! /MARGIN_RIGHT=-
! /PAGE_NUMBER_POSITION=-
! /PAGE_NUMBER_ALIGNMENT=-
! /PAGE_TAG_LINE_2=-
! /PAGE_TAG_LINE_3=-
! /PAGE_TAG_LINE_4=-
;
```

Manually Set Report Options

These parameters can be changed in the pipeline command, or changed manually by right clicking on the report page and selecting **Report Options**.

Page Layout

Page Layout Tab:

/PAGE_LAYOUT= Portrait or Landscape

/PAPER_SIZE= Letter, Legal or A4

Lab Information

The Lab Information can be added in the Laboratory Information tab. Laboratory Information Tab: /COLLECTOR= /FACILITY_LINE_2= /ADDRESS_LINE_1= /ADDRESS_LINE_2= /ADDRESS_LINE_3= /PHONE= /FAX= /NOTES_LINE_1= /NOTES_LINE_2= /NOTES_LINE_3= /LOGO=	ReportOptions.jpeg { {/images/thumb/a/a7/ReportOptions.jpeg/500px-ReportOptions.jpeg
Note that to see the lab information in the report, the lab information must be “graphed”. By default facility line 1 is set to the laboratory name and this parameter cannot be modified.	ReportOptions_LabInfo_Report.jpeg { {/images/thumb/c/c3/ReportOptions_LabInfo_Report.jpeg/500px-ReportOptions_LabInfo_Report.jpeg

Page Layout

Page Layout Tab: 1. Margins: /MARGIN_TOP= /MARGIN_BOTTOM= /MARGIN_LEFT= /MARGIN_RIGHT= 2. Page Number Format: /PAGE_NUMBER_POSITION= /PAGE_NUMBER_ALIGNMENT= Note paper size and orientation can be set in the first two parameters of this pipeline command.	ReportOptions_PageLayout.jpg { {/images/thumb/e/ef/ReportOptions_PageLayout.jpg/500px-ReportOptions_PageLayout.jpg
---	---

Printing

Printing Tab: /PAGE_TAG_LINE_2= /PAGE_TAG_LINE_3= /PAGE_TAG_LINE_4=	ReportOptions_Printing.jpg { {/images/thumb/b/b7/ReportOptions_Printing.jpg/500px-ReportOptions_Printing.jpg
---	---

From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:report_commands:set_report_options&rev=1718648556

Last update: 2024/06/17 18:22

