

the lab information can be added in the laboratory information tab.\\ \\ \\ **laboratory information tab:**\\ \\ **/collector=**\\ **/facility_line_2=**\\ **/address_line_1=**\\ **/address_line_2=**\\ **/address_line_3=**\\ **/phone=**\\ **/fax=**\\ **/notes_line_1=**\\ **/notes_line_2=**\\ **/notes_line_3=**\\ **/logo=**\\ \\ | reportoptions.jpeg\\ | |\\ \\ note that to see the **lab information** in the report, the lab information must be "graphed".\\ \\ by default facility line 1 is set to the laboratory name and this parameter cannot be modified.\\ \\ | reportoptions_labinfo_report.jpeg\\ | ===== page layout ===== |**page layout tab:**\\ \\ 1. margins:\\ **/margin_top=**\\ **/margin_bottom=**\\ **/margin_left=**\\ **/margin_right=**\\ \\ 2. page number format:\\ **/page_number_position=**\\ **/page_number_alignment=**\\ \\ note paper size and orientation can be set in the first two parameters of this pipeline command.\\ | reportoptions_pagelayout.jpg| ===== printing ===== |**printing tab:**\\ \\ **/page_tag_line_2=**\\ **/page_tag_line_3=**\\ **/page_tag_line_4=**\\ | reportoptions_printing.jpg| } } } } }

From:

<https://has-motion.com/wiki/> - **Wiki Documentation Site**

Permanent link:

https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:report_commands:set_report_options&rev=1718804324

Last update: **2024/06/19 13:38**

