

Set Report Page title

Set the title of a report page.

Edit Report Page Title

×

Enter the title that you would like on page:

Arial

▼

12

▼

Black

▼

Custom

B

I

U

≡

≡

≡

☐ Evaluate title as an expression before display

OK

Cancel

```
Set_Report_Page_Title
! /PAGE_NUMBER=
! /PAGE_TITLE=
;
```

From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:report_commands:set_report_page_title

Last update: 2024/07/17 15:46

