

Set the title of a report page.



```
Set_Report_Page_Title
! /PAGE_NUMBER=
! /PAGE_TITLE=
;
```

Retrieved from ""

From:  
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:  
[https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:report\\_commands:set\\_report\\_page\\_title&rev=1718644097](https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:report_commands:set_report_page_title&rev=1718644097)

Last update: 2024/06/17 17:08

