

Delete Folder

The Delete_Folder command allows users to delete specific folders of signals or metrics from the [data tree](#).

Pipeline Command

The Delete_Folder command is located in the pipeline workshop's Signal/Data Management folder.

```
Delete_Folder
/SIGNAL_TYPES=
/FOLDER=
! /REMOVE_FROM_RECALC_PIPELINE=FALSE
;
```

Command Parameters

Parameter	Description
Signal Types	The signal type to work with.
Folder	The specific folder(s) to be deleted.
Remove from Recalc pipeline	If true, then the signals contained in the deleted folder will also be removed from the Recalc pipeline . Default value is false.

Example: Remove all Link_Model_Based items

This example demonstrates how to remove all [LINK_MODEL_BASED](#) items from the workspace and the Recalc pipeline.

```
Delete_Folder
/SIGNAL_TYPES=LINK_MODEL_BASED+LINK_MODEL_BASED
/FOLDER=ORIGINAL+PROCESSED
/REMOVE_FROM_RECALC_PIPELINE=TRUE
;
```

After running this command, all ORIGINAL and PROCESSED Link_Model_Based signals will be removed from the workspace and the corresponding [Compute_Model_Based_Data](#) commands will be removed from the Recalc pipeline.

From:
<https://has-motion.com/wiki/> - **Wiki Documentation Site**

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:signal_commands:delete_folder

Last update: **2026/09/09 18:43**

