

remove a signal from the data tree.

the details on the command are found below:

```
\\ remove_signals\\ /signal_types=the type of signal to be evaluated\\  
/signal_folder==the signal folder\\ /include_models= true or false assigning  
true will remove signals from models assigned to active files\\
```

examples

for example, remove one derived signal and one metric

remove_signals /signal_types=derived+metric /signal_names= derived1+metric1 !
/signal_folder=processed+processed ; signals can also be removed manually from the data tree

[data_tree#deleting_a_signal](#)

From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:signal_commands:remove_signals&rev=1718801598

Last update: 2024/06/19 12:53

