

2026/09/16 23:22	1/1	Second Derivative of a Signal
Second_Derivative		
/Signal_Types=	The type of signal to be processed	
/Signal_Folder=	The name of the signal folder	

Compute the second derivative at each point of a signal. The first and last data points become DATA_NOT_FOUND(identified in the C3D file as a POINT Reliability value of &endash;1).

Derivatives are calculated using Finite Difference Algorithms.

Given the signal:

x(ti) for i=1, 2, ..., n The second derivative is calculated using: $xDDot(ti)= (x(ti+1)-2x(ti)+x(ti-1)) / ((ti+1)-(ti))^2$

From: <https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link: https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:signal_commands:second_derivative_of_a_signal&rev=1720791093

Last update: 2024/07/12 13:31

