

✖ 2026/09/17 04:07 Second_Derivative	1/1	Second Derivative of a Signal
/Signal_Types=	The type of signal to be processed	
/Signal_Folder=	The name of the signal folder	

Compute the second derivative at each point of a signal. The first and last data points become DATA_NOT_FOUND(identified in the C3D file as a POINT Reliability value of &endash;1).

Derivatives are calculated using Finite Difference Algorithms.

Given the signal:

$x(t_i)$ for $i=1, 2, \dots, n$ The second derivative is calculated using: $xDDot(t_i)= (x(t_{i+1})-2x(t_i)+x(t_{i-1}))/((t_{i+1})-(t_i))^2$

From:

<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:

https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:signal_commands:second_derivative_of_a_signal&rev=1720793235

Last update: 2024/07/12 14:07

