

| \ the

[[visual3d:documentation:pipeline:report_commands:add_graph_annotation|add_graph_annotation]]
command specifies which graph to annotate using the draw_order parameter. \ \ ****the order the
graphs were drawn in can be determined in two ways:**** \ \ 1) the order the graphs are listed on the
legend (shown to the right, highlighted in red). the first graph listed has a draw_order of 1 \ \ 2) the
order the graphs are listed in the report editor (shown below highlighted in red). the first graph listed
has a draw_order of 1 \ \ the order the graphs are drawn can be changed by using the arrows to the
right of the scroll bar in the list of report items and saving the new report template. (below) \ \
graphannotations_draworder_reportlist.jpg |graphannotations_draworder_legend.jpg | }

From:

<https://has-motion.com/wiki/> - **Wiki Documentation Site**

Permanent link:

https://has-motion.com/wiki/doku.php?id=visual3d:documentation:reports:draw_order&rev=1718801613

Last update: **2024/06/19 12:53**

