

Draw Order

The orderThe graphs were drawn in can be determined in two ways:

- 1) The order the graphs are listed on the legend (shown to the right, highlighted in red). The first graph listed has a Draw_Order of 1
- 2) The order the graphs are listed in the report editor (shown below highlighted in red). The first graph listed has a Draw_Order of 1

The order the graphs are drawn can be changed by using the arrows to the right of the scroll bar in the list of report items and saving the new report template. (below)



From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:reports:draw_order&rev=1720027819

Last update: 2024/07/03 17:30

