

# Draw Order

The order the graphs were drawn in can be determined in two ways:

- 1) The order the graphs are listed on the legend (shown to the right, highlighted in red). The first graph listed has a Draw\_Order of 1
- 2) The order the graphs are listed in the report editor (shown below highlighted in red). The first graph listed has a Draw\_Order of 1

The order the graphs are drawn can be changed by using the arrows to the right of the scroll bar in the list of report items and saving the new report template. (below)



From:  
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:  
[https://has-motion.com/wiki/doku.php?id=visual3d:documentation:reports:draw\\_order&rev=1721149239](https://has-motion.com/wiki/doku.php?id=visual3d:documentation:reports:draw_order&rev=1721149239)

Last update: 2024/07/16 17:00

