

Draw Order

The order the graphs were drawn in can be determined in two ways:

- 1) The order the graphs are listed on the legend (shown to the right, highlighted in red). The first graph listed has a Draw_Order of 1
- 2) The order the graphs are listed in the report editor (shown below highlighted in red). The first graph listed has a Draw_Order of 1

The order the graphs are drawn can be changed by using the arrows to the right of the scroll bar in the list of report items and saving the new report template. (below)

Workspace | Signals and Events | Models | Reports | CalTester | Real-Time


Pg	Col	Row	Type	Name	Subtype	Comp	T
1	1	1	TEMPORAL_DISTA...				A
2	1	1	LINK_MODEL_BAS...	Left_Knee_Angle	JOINT_ANGLE	X	C
2	1	1	LINK_MODEL_BAS...	Left_Knee_Angle	JOINT_ANGLE	X	C

Modify Selected | Export Selected | Delete Selected

Add a new report item to the report

Page 2 | Column 1 | Span 1 | Item to add: 2D Graph | Add

Row 1 | Span 1



From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:reports:draw_order&rev=1721229096

Last update: 2024/07/17 15:11

