

Language: English • français • italiano • português • español ****

==== Contents ====

* 1 Introduction

* 2 To change the range:

Introduction

How to Modify Explicit Settings for Report Graphs

By default, the Y axis of a graph will be graphed from max to min:

Changing the settings will allow you to manually specify the max/min value:

Consider a graph created with the settings on the right:
Since the Y-Axis Scale is set to Global Min to Max, the graph that is created has the range of data for the y-axis determined by the maximum and minimum values of the signals graphed:

To change the range:

1. In the Reports tab, select the graph using the Left Mouse Button

- Handles will appear indicating that the graph is active

- If the handles don't appear, the graph is not selected!!

2. Right click on the selected graph

3. Select Format Graph

4. Set the y-axis range:

- Check Modify Y Axis Scale

- Specify the Y Min/Y Max

5. Click Apply

The graph has now been explicitly set to these properties:

The Report Graph dialog will now reflect that the settings are explicit. If the graph template is saved, these explicit settings will be saved in the **Report Template**



Retrieved from ""

From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:reports:explicit_settings_for_report_graphs&rev=1718386247

Last update: 2024/06/14 17:30

