

Explicit_Settings_for_Report_Graphs

Introduction

How to Modify Explicit Settings for Report Graphs

By default, the Y axis of a graph will be graphed from max to min: 	Changing the settings will allow you to manually specify the max/min value: 
---	--

Consider a graph created with the settings on the right: Since the Y-Axis Scale is set to Global Min to Max, the graph that is created has the range of data for the y-axis determined by the maximum and minimum values of the signals graphed: 	
--	--

To change the range:

1. In the Reports tab, select the graph using the Left Mouse Button - Handles will appear indicating that the graph is active - If the handles don't appear, the graph is not selected!! 2. Right click on the selected graph 3. Select Format Graph	
---	--

4. Set the y-axis range: - Check Modify Y Axis Scale - Specify the Y Min/Y Max 5. Click Apply	
---	--

The graph has now been explicitly set to these properties:



The Report Graph dialog will now reflect that the settings are explicit. If the graph template is saved, these explicit settings will be saved in the Report Template	
--	--

From:
<https://has-motion.com/wiki/> - **Wiki Documentation Site**

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:reports:explicit_settings_for_report_graphs&rev=1720793252

Last update: **2024/07/12 14:07**

