

add_graph_annotation]] annotationdlg1.jpg the initial dialog. select **add** to add an annotation to the graph select **edit** to modify an existing annotation. annotationdlg2.jpg the initial dialog for the annotation annotationdlg3.jpg select the parameters of the annotation. if all graphs are selected an annotation is added for every occurrence of the event for every graph. annotationgraph1.jpg notice the multiple vertical lines because there is an event for each graph. annotationdlg4.jpg if average event or numerical value is selected, visual3d interprets the edit box. for example, specifying the event label t3 annotationgraph2.jpg for example, in this case the average location of the event is identified and drawn. annotationdlg5.jpg if average event or numerical value is selected, visual3d interprets the edit box. for example, specifying the time 3 annotationgraph3.jpg for example, in this case the annotation is drawn at 25% of the x-axis..

==== section 3 ==== the signal used for the x-axis. reportgraph4.jpg\\ in this example, time is displayed on the x-axis. the signals in the y-axis properties for (file to graph and within a graph) that have occurrences of the range range begin to **range end** will be drawn on the graph. if **intermediate events** are selected, only ranges from **range begin** to **range end** that contain the **intermediate events** will be graphed. if exclude sequences containing is populated, only ranges from **range begin** to **range end** that do not contain the **excluded events** will be graphed

==== section 4 ==== the signal used for the y-axis reportgraph5.jpg in this example, the x-component of the right ankle angle will be displayed. the range of the graph can be specified as:

- normalize 0% to 100% - each range of data will be scaled from 0 to 100 percent
- global min to max - the minimum and maximum value of the signals across all ranges is computed
- by explicit setting - the user can format the graph and specify an explicit range regardless of the signals selected.

=== defining an explicit range === select a graph in the report. if the graph is selected you will see "handles" reportgraph6.jpg move the cursor over the selected graph, and push the right mouse button. a menu will appear. reportgraph7.jpg select format graph and a dialog will appear reportgraph8.jpg select the y min and y max values and save these settings as the default. this means that these setting will be preserved in the report template. reportgraph9.jpg

==== section 5 ====

- === set the properties of the legend === move the cursor over the selected graph, and push the right mouse button. a menu will appear. reportgraph7.jpg select format graph and a dialog will appear reportgraph8.jpg select the option to display the legend. select the button to modify the legend reportgraph12.jpg the legend can be placed at the top left, top right, bottom left, or bottom right. 1 or more of the selections of the legend text can be selected.
- === adding custom legend text === adding the text for the report legend. reportgraph9a.jpg if you have populated this dialog, you will not see the legend in the graph by default. if you have entered text in the graph format dialog, you need to select the option to display text in the "legend properties" dialog.
- ===== adding more signals to the report graph ===== the report dialog acts on one report graph entity. the graph that appears in the report, however, may contain many entities. a second signal can be added to the graph in the report by adding another graph entity to the same location. each graph entity may contain its own graph properties. reportgraph10.jpg\\ for example, this image shows two graphs rankle_angle and lankle_angle located on page 2, column 1, row 1. the resulting graph contains data from the 2 signals. rangle_angle in red and lankle_angle in black, each with its own mean and standard deviation. reportgraph11.jpg

}}

From: <https://has-motion.com/wiki/> - **Wiki Documentation Site**

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:reports:report_2d_graph&rev=1718804354

Last update: **2024/06/19 13:39**

