

Language: English • français • italiano • português • español ****

==== Contents =====

- * 1 Table Formatting
 - * 1.1 Metric Table Size
 - * 1.2 Cell Borders
 - * 1.3 Select Area
 - * 1.4 Metric Fonts
- * 2 Metric Table Items
 - * 2.1 Metric Data item
 - * 2.2 Subject Data Item
 - * 2.3 Text Item

Text Metric Table



Table Formatting

Metric Table Size



Cell Borders



Select Area



Metric Fonts



Metric Table Items

Metric Data item



The combo box will contain a list of all text signals and metric signals in the open Workspace

For example:

Create a text item containing the c3d parameter POINT UNITS

[Set_Pipeline_Parameter_To_Data_Value](#) /SIGNAL_TYPES=PARAMETERS /SIGNAL_FOLDER=POINT /SIGNAL_NAMES=UNITS ! /SIGNAL_COMPONENTS=ALL_COMPONENTS /PARAMETER_NAME=TEXT ;

Create a text signal using the pipeline parameter

[Create_Text_Data](#) /SIGNAL_NAME=SCOTT /SIGNAL_FOLDER=TEST /SIGNAL_DATA_VALUE=::TEXT ;

This Text signal is available in the Table.



Subject Data Item



The edit box will contain all subject data from the current Workspace

For example, to add the Subject's Name to the metric table, you can define a Subject Metric to be a [text string](#).

Text Item

Enter the text directly into the edit box.



Retrieved from ""

From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:reports:text_metric_table&rev=1718386259

Last update: **2024/06/14 17:30**

