

Text Metric Table



Table Formatting

Metric Table Size



Cell Borders



Select Area



Metric Fonts



Metric Table Items

Metric Data item



The combo box will contain a list of all text signals and metric signals in the open Workspace

For example:

Create a text item containing the c3d parameter POINT UNITS

[Set Pipeline Parameter To Data Value](#) /SIGNAL_TYPES=PARAMETERS /SIGNAL_FOLDER=POINT /SIGNAL_NAMES=UNITS ! /SIGNAL_COMPONENTS=ALL_COMPONENTS /PARAMETER_NAME=TEXT ;
Create a text signal using the pipeline parameter

[Create Text Data](#) /SIGNAL_NAME=SCOTT /SIGNAL_FOLDER=TEST /SIGNAL_DATA_VALUE=::TEXT ;
This Text signal is available in the Table.



Subject Data Item



The edit box will contain all subject data from the current Workspace

For example, to add the Subject's Name to the metric table, you can define a Subject Metric to be a [text string](#).

Text Item

Enter the text directly into the edit box.



From:
<https://has-motion.com/wiki/> - **Wiki Documentation Site**

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:reports:text_metric_table&rev=1718717673

Last update: **2024/06/18 13:34**

