

Matlab_FAQ

This page is the FAQs regarding Matlab.

FAQs

Question: I am having trouble importing data compiled in Matlab into Visual 3D. What is the most likely problem?

Answer (click to expand)

When exporting from matlab, you need to use the option to save an old format (uncompressed) .mat file. The [Matlab](#) page has specific instructions for importing and exporting from Matlab.

From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:third-party:matlab:matlab_faq&rev=1720793268

Last update: **2024/07/12 14:07**

