

Natural_Point_Overview

[Natural Point\]](#)

[OptiTrack\]](#)

The developers of the OptiTrack cameras. These cameras use their Motive software to capture data and save C3D files, and can stream real-time data to Visual3D.

The camera capture software included with [AMASS](#) uses the NaturalPoint SDK to capture data, and AMASS produces C3D files.

From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:third-party:natural_point:natural_point_overview&rev=1721149356

Last update: **2024/07/16 17:02**

