

Natural Point Overview

Natural Point are the developers of the OptiTrack cameras.

These cameras use their Motive software to capture data and save C3D files, and can stream real-time data to Visual3D. The camera capture software included with AMASS uses the NaturalPoint SDK to capture data, and AMASS produces C3D files.

From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:third-party:natural_point:natural_point_overview&rev=1731690485

Last update: 2024/11/15 17:08

