

Natural Point Overview

Natural Point are the developers of the **OptiTrack** cameras.

These cameras use their Motive software to capture data and save C3D files, and can stream real-time data to Visual3D. The camera capture software included with **AMASS** uses the NaturalPoint SDK to capture data, and AMASS produces C3D files.

Plugins

Has-Motion offers multiple plugins for real-time data streaming into Visual3D. When using a real-time streaming plug-in with Optitrack camera systems, the type of camera does not change how the data is streamed into Visual3D. Optitrack processes their camera data the same regardless of version before sending data to Visual3D. If timestamps are needed for Optitrack data while using a real-time streaming plugin, an **esync2** will be needed.

From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:third-party:natural_point:natural_point_overview&rev=1757611062

Last update: 2025/09/11 17:17

