




An example:

Workspace History

Open

Date: 20230721 Time: 102721
Visual3D Version: 1.0000
CMD: C:\Users\amyco\Documents\Inspect3D Regression Test (2023-07-19)\Inspect3D Regression Test\Data_Original\Sub

Save

Date: 20230721 Time: 102734
Visual3D Version: 1.0000
CMD: C:\Users\amyco\Documents\Inspect3D Regression Test (2023-07-19)\Inspect3D Regression Test\Data_Original\Sub

Open

Date: 20230721 Time: 102806
Visual3D Version: 1.0000
CMD: C:\Users\amyco\Documents\Inspect3D Regression Test (2023-07-19)\Inspect3D Regression Test\Data_Original\Sub

Save

Date: 20230721 Time: 102818
Visual3D Version: 1.0000
CMD: C:\Users\amyco\Documents\Inspect3D Regression Test (2023-07-19)\Inspect3D Regression Test\Data_Original\Sub

Open

Date: 20250506 Time: 143012
Visual3D Version: 1.0001
CMD: C:\Users\amyco\Documents\Data_Original\Sub02\Sub02_Workspace.cmz

Save

Date: 20250506 Time: 143013
Visual3D Version: 1.0001
CMD: C:\Users\amyco\Documents\Data_Original\Sub02\Sub02_Workspace.cmz

Open

Date: 20250528 Time: 155328
Visual3D Version: 2025.050200
CMD: C:\Users\amyco\Documents\Data_Original\Sub02\Sub02_Workspace.cmz

OK

From:
<https://has-motion.com/wiki/> - **Wiki Documentation Site**

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:getting_started:workspace_history&rev=1721230998

Last update: **2024/07/17 15:43**

