


Real Time Interface

Introduction to the Visual3D RealTime Interface Tab



Visual3D Real Time interface tab
as of March 2022

- Stream Information:
 - Built-in RT vs External RT
 - Camera Frame Rate
 - Built-In RT Support
 - Use FP Parameters from CalFile
 - Force Platform Zero Baseline
 - Hardware Settings
- Real Time Pipeline Commands
 - List of available pipeline commands that run in real time
- Real-Time Graphs

- 2D/3D Graph Toggle
- Graph Display Plane
- Full Screen Display
- Play Sounds
- Switch Vertical Axis (+/-)
- Switch Horizontal Axis (+/-)
- Frames to Graph
 - Determines how many frames the real-time keeps in the data buffer
 - If you want to plot 10 gait cycles, the number of frames must include enough frames to plot 10 gait cycles
- Max Frame Gap
- Define graphs

From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:real_time_streaming:real_time_interface&rev=1721229830

Last update: **2024/07/17 15:23**

