

Language: **English** • [français](#) • [italiano](#) • [português](#) • [español](#) ****

===== Contents =====

- * [1 Landmark Type 1](#)
- * [2 Landmark Type 2](#)
- * [3 Landmark Type 3](#)
- * [4 Landmark Type 4](#)
- * [5 Landmark Type 5](#)
- * [6 Landmark Type 6](#)
- * [7 Landmark Type 7](#)
- * [8 Landmark Type 8](#)
- * [9 Automatically Created Hip Landmarks](#)
- * [10 Functional Joint Landmarks](#)
- * [11 Tutorial: Digitizing Pointer Post Processing](#)

This Tutorial page lists the many different ways to create Landmarks within Visual3D.

Landmark Type 1

Creating a landmark at the end of a segment that is tracked with the segment.

Landmark Type 2

Projecting a marker onto a plane.

Landmark Type 3

Projecting a marker onto one of the principal planes of the laboratory.

Landmark Type 4

Projecting a marker onto a line.

Landmark Type 5

Creating a landmark that defines the location of a segment endpoint relative to 3 tracking markers (or landmarks).

Landmark Type 6

Create a landmark relative to a Segment Coordinate System.

Landmark Type 7

Create a landmark relative to the Global Coordinate System (Lab).

Landmark Type 8

Automatically Created Hip Landmarks

Landmarks defining the location of the hip joint center automatically created when the CODA or Helen Hayes Pelvis segment is created.

Functional Joint Landmarks

Calculation of a functional joint requires movement of one segment relative to another segment. The algorithm searches for a point (or for a one degree of freedom joint, an axis) that is stationary relative to the 2 segments (or 2 sets of markers).

Tutorial: Digitizing Pointer Post Processing

Identifying anatomical landmarks without having to place motion capture markers at the location using a digitizing pointer.

Retrieved from ""

From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:tutorials:modeling:creating_landmarks&rev=1718386346

Last update: 2024/06/14 17:32

