

This Tutorial page lists the many different ways to create Landmarks within Visual3D.

## **Landmark Type 1**

Creating a landmark at the end of a segment that is tracked with the segment.

## **Landmark Type 2**

Projecting a marker onto a plane.

## **Landmark Type 3**

Projecting a marker onto one of the principal planes of the laboratory.

## **Landmark Type 4**

Projecting a marker onto a line.

## **Landmark Type 5**

Creating a landmark that defines the location of a segment endpoint relative to 3 tracking markers (or landmarks).

## **Landmark Type 6**

Create a landmark relative to a Segment Coordinate System.

## **Landmark Type 7**

Create a landmark relative to the Global Coordinate System (Lab).

## **Landmark Type 8**

## Automatically Created Hip Landmarks

Landmarks defining the location of the hip joint center automatically created when the CODA or Helen Hayes Pelvis segment is created.

## Functional Joint Landmarks

Calculation of a functional joint requires movement of one segment relative to another segment. The algorithm searches for a point (or for a one degree of freedom joint, an axis) that is stationary relative to the 2 segments (or 2 sets of markers).

## Tutorial: Digitizing Pointer Post Processing

Identifying anatomical landmarks without having to place motion capture markers at the location using a digitizing pointer.

From:  
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:  
[https://has-motion.com/wiki/doku.php?id=visual3d:tutorials:modeling:creating\\_landmarks&rev=1719432500](https://has-motion.com/wiki/doku.php?id=visual3d:tutorials:modeling:creating_landmarks&rev=1719432500)

Last update: **2024/06/26 20:08**

