

# Creating\_Landmarks

This Tutorial page lists the many different ways to create Landmarks within Visual3D.

## Landmark Type 1

Creating a landmark at the end of a segment that is tracked with the segment.

## Landmark Type 2

Projecting a marker onto a plane.

## Landmark Type 3

Projecting a marker onto one of the principal planes of the laboratory.

## Landmark Type 4

Projecting a marker onto a line.

## Landmark Type 5

Creating a landmark that defines the location of a segment endpoint relative to 3 tracking markers (or landmarks).

## Landmark Type 6

Create a landmark relative to a Segment Coordinate System.

## Landmark Type 7

Create a landmark relative to the Global Coordinate System (Lab).

## Landmark Type 8

## Automatically Created Hip Landmarks

Landmarks defining the location of the hip joint center automatically created when the CODA or Helen Hayes Pelvis segment is created.

## Functional Joint Landmarks

Calculation of a functional joint requires movement of one segment relative to another segment. The algorithm searches for a point (or for a one degree of freedom joint, an axis) that is stationary relative to the 2 segments (or 2 sets of markers).

## Tutorial: Digitizing Pointer Post Processing

Identifying anatomical landmarks without having to place motion capture markers at the location using a digitizing pointer.

From:  
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:  
[https://has-motion.com/wiki/doku.php?id=visual3d:tutorials:modeling:creating\\_landmarks&rev=1720027845](https://has-motion.com/wiki/doku.php?id=visual3d:tutorials:modeling:creating_landmarks&rev=1720027845)

Last update: 2024/07/03 17:30

